



## Portable Freshness

### Oral Appliance Device Cleaner

Colour: Light Blue

Scent: Mint

#### Available Sizes:

60 mL Travel Sprayer (case of 10) (05-MOFF-060)

#### Purpose:

MouthOFF™ cleans and refreshes oral appliances without staining, while minimizing the impact on the Environment, People and Equipment.

#### Suitable for:

Cleaning and protecting oral appliances / devices including dentures, sleep apnea machines, mouthguards, soft plastics, metals, silicone and acrylics. Compatible with Nebulizers, Aerochamber® and other devices with silicone masks / mouthpiece.

#### Applied Cleaning Concepts:

MouthOFF™ contains biodegradable surfactants for cleaning and protecting oral appliances.

#### Quality Commitment:

Our vision of cleaning holistically links the patient/client safety with the safety of professional staff. Our chemistry is based on materials of pharma-grade (USP EP BP NF) origin. Quality Ingredients include biodegradable surfactants, chlorhexidine, xylitol with USP grain alcohol and high purity water. OECD 301D Tests conclude biodegradability (28 days).

#### Instructions for Use on Oral Appliances:

Rinse the device thoroughly. Spray or submerge the device for a minimum of 5 min. Rinse the device again before use.

#### Instructions for Use on Sleep Apnea Machines:

- 1) Unplug your CPAP machine from the power source.
- 2) Disassemble your mask into 3 parts (headgear, cushion and frame).
- 3) Cushion should be cleaned daily using MouthOFF™. Spray MouthOFF™ generously on the cushion and rinse thoroughly. Dry cushion with a clean Micrylium LeCLOTH™. Headgear and frame should be cleaned weekly with Micrylium BioSON™. In a sink or tub, clean your mask cushion and headgear to remove any oils. Gently rub with a mild liquid detergent and warm, drinking-quality water. Avoid using stronger cleaning products as they may damage the mask or leave harmful residue.
- 4) Use BioSON™ at a dilution rate 1:200 which is 5mL of solution for 1L of potable water in a sink or tub to clean and remove oils from the mask cushion and headgear. Gently scrub the appliance in the solution with a saturated LeCLOTH™. Simultaneously immerse the air tubing in the solution.
- 5) Rinse all components thoroughly with warm potable water being sure to remove all residue.
- 6) Place on a flat surface, on top of a towel, to dry. Use a dry LeCLOTH™ to remove any moisture from irregular areas where it may pool. Avoid placing in direct sunlight.

#### Precautions:

USE FULL STRENGTH – Do not dilute. RTU (Ready To Use Formulation)

Contains 9.5% Ethanol and Xylitol. (Keep away from Pets)

#### Helpful Tips:

To ensure good results, rinse the device thoroughly before application.

See our website for specific details: [www.micrylium.com/products](http://www.micrylium.com/products)



Revision Date: Dec 15, 2023 v2.0




KINDER



SAFER



FASTER



1. IDENTIFICATION									
Product Name		Mouth Off			Manufacturer		Micrylium Laboratories Inc.		
Registration		EU CE	Class I		Address		5000M Dufferin Street, Toronto, Canada, M3H 5T5 www.micrylium.com		
		EU CE BAG	Class I						
Indication		Medical device cleaner			Phone		416-667-7040		
Emergency Phone #		CHEMTREC			Fax		416-667-0071		
					1-800-424-9300		CANUTEC		1-613-996-6666
2. HAZARD IDENTIFICATION									
Symbol Pictogram		Not Applicable			Signal Word		Not Applicable		
					Symbol		Not Applicable		
Classification		Not Applicable							
Health Hazard		No Serious Health Hazards				Environmental Hazards		Biodegradeable	
Precautionary & Hazard Statements		P102: Keep out of reach of children. P301: IF SWALLOWED: Drink large quantities of water or milk. P305: IF IN EYES: Flush eyes with large quantities of water.			H302: Harmful if swallowed.				
3. COMPOSITION									
Chemical		CAS #		LD-50 (Oral, mg/kg)(Rat)		Concentration (%)			
Ethanol		64-17-5		7,060		9.5%			
Chlorhexidine Gluconate		18472-51-0		2,000		0.12%			
Xylitol		87-99-0		4,000		0.6%			
4. FIRST AID MEASURES									
Inhalation		If breathing is difficult, remove individual to fresh air.			Ingestion		Drink large quantities of milk or water to dilute. Do not induce vomiting.		
Skin Contact		No adverse effects.			Eye Contact		Flush with plenty of water.		
Most Important Symptoms and Effects (Acute and Delayed)									
May cause mild reversible irritation if there is eye contact or is inhaled.									
Indication of any Immediate Medical Attention and Special Treatment Needed									
Not Applicable.									
5. FIREFIGHTING MEASURES									
Non-flammable.									
6. ACCIDENTAL RELEASE									
Use all means to prevent spillage.									
7. HANDLING & STORAGE									
Store in a cool, dry, well ventilated location. Keep away from heat, sparks and flames. DO NOT mix with bleach or peroxides. Storage and Transport: 5°- 30°C									
8. EXPOSURE CONTROLS/ PERSONAL PROTECTION									
No specific measures required.									
9. PHYSICAL AND CHEMICAL PROPERTIES									
Physical State	Colour	Odour	Solidification point	Boiling point	Flash Point	Density	pH	Kinematic Viscosity	
Transparent, Liquid	Light Blue	Peppermint	-5°C	105°C	N/A	6.5	9.0	12 mm²/s	
10. STABILITY AND REACTIVITY									
Stable under normal conditions. Incompatibility: Strong oxidants, acid chlorides, silver salts Decomposition: Products: CO <sub>2</sub> , CO									
11. TOXICOLOGICAL DATA									
Acute Dermal Toxicity		Not found to be dermal sensitizer			Acute Oral		N/A		
Ocular Irritation		N/A			Acute Inhalation Toxicity		N/A		
Reproductive Hazards		Ingestion/inhalation can be harmful. (TDLo 300mg/Kg Ethanol)			Carcinogenicity		Ingestion of Ethanol IARC Group 1.		
CAUTION: Xylitol is toxic to dogs, rabbits cows, baboons and horses when consumed.									
12. ECOLOGICAL INFORMATION									
Surfactants are readily biodegradable in soil and water. Persistence unlikely based on available data.									
Ethanol		EC50 (72h) = 275 mg/l (Chlorella vulgaris)		Fathead minnow (Pimephales promelas) LC50 = 14200 Mg/L/96h		Photobacterium Phosphoreum:EC50 = 34634 Mg/L/30 min Photobacterium Phosphoreum:EC50 = 35470 Mg/L/5 min		EC50 = 9268 mg/L/48h EC50 = 10800 mg/L/24h	
13. DISPOSAL CONSIDERATIONS									
Domestic.									
14. TRANSPORT INFORMATION									
Not regulated									
N/A	Land			Sea			Air (IATA)		
	N/A			N/A			N/A		
	N/A			N/A			N/A		
	N/A			N/A			N/A		
15. REGULATORY INFORMATION									
TSCA – No reporting required.					CERCLA – No hazardous pollutants or ozone depletion.				
16. OTHER INFORMATION									
The information and recommendations contained herein are based on information believed to be correct. It is offered in good faith, without guarantee. Micrylium Laboratories Inc. make no warranty expressed or implied.									
Effective Date: 2022/10/27			Revision Date: 2025/10/27			Document: MouthOFF 1.3			